

Micro-nanofabricated platform technology for cell seeding experiments in neuro-nanobiology

R. Luttge

*MESA+ Institute for Nanotechnology, University of Twente, 7500AE Enschede,
The Netherlands
r.luttge@utwente.nl*

Pioneering research in the neurosciences utilizes micro-nanotechnology.¹ The exact definition of cell seeding environments and tuning of scaffold properties for dissociated primary neurons, however, remains a challenge. By means of a polydimethylsiloxane (PDMS) soft lithography, a micro-nanofabricated silicon template is applied to provide for a new cell culture platform technology, here introduced for experiments in neuro-nanobiology. The PDMS culture platforms contain micro-nanoscale structured wells offering cell guidance behavior. The use of similar silicon templates for nanobiology has been demonstrated earlier in the study of bone tissue regenerative medicine.² As in the previous work, a Laser Interference Lithography (LIL) pattern transfer process was applied to produce a nanostructure across the surface of a 4" prime silicon wafer and subsequent PDMS replication resulted into multiple micro-nanostructured cell seeding platforms. Figure 1 A and B show the micro/nanostructured silicon template and Figure 1 C, D, and E depict some examples of the PDMS replicated structures after cell seeding and network growth. Micro-nanostructured PDMS surfaces were assembled with rings also from PDMS to provide a larger cell seeding reservoir to contain the initial culture medium. To carry out a pilot cell culture repetitive micro-nanostructured platforms were inserted into the 24 wells of a standard tissue culture plate and primary neurons dissociated from the brains of newborn rats were seeded and visually followed over 8 weeks. In-vitro cell culture medium was refreshed every two days. Figure 1 D and E specifically show that the neurons are sustainable during prolonged culture on the PDMS micro-nanostructured platforms. The network-forming growth behavior of the primary neurons were preserved after seeding as staining the cells with a LIVE/DEAD® viability/cytotoxicity kit for mammalian cells after completion of 8 weeks confirmed that a majority of cells were flourishing.³

Thanks are due to M. Domanski, S. Schlautmann, E. Landman and C.J.M. Groot Jebbink-van Ginkel for assistance in fabrication and neuronal cell handling. Laboratory access to University of Twente's MIRA Institute for Biomedical Engineering and Technical Medicine and financial support by STW are greatly acknowledged.

¹ Alberti M., et al. (2010) *Microfluid Nanofluid* 9:963-972

² Lamers, E., et al. (2010) *Biomaterials*, 31 (12), pp. 3307-3316.

³ <http://tools.invitrogen.com/content/sfs/manuals/mp03224.pdf>




Figure 1: Micro-nanostructured silicon template for replication of PDMS cell seeding platforms for neuro-nanobiology: A. Pillar of microscale dimensions with a diameter of $600\mu\text{m}$ and a height of $100\mu\text{m}$ and a nanostructured top forming subsequently a nanostructured well in the PDMS cell seeding platform. B. Detail of the pillar top surface carrying a well-defined nanograting structure with a periodicity of 300 nm . C. Optical microscope image of a replicated PDMS platform with primary neurons seeded onto its surfaces. On both surfaces, atop the well-structure and inside of the wells neurons form clustered networks. The micro-nanostructured wells form a barrier to excessive cellular network growth, while the unstructured platform surface provides carpet-like cellular network growth. The neuronal clusters inside of the well are interconnected with the neurons atop. D. Fluorescent optical image of two wells after 8 weeks of cell culturing depicting neuronal clusters after a live/dead staining. E. Detail of fluorescent image inside of one of the wells.