

Enhanced Removal of Material by FIB Etching at Glancing Angle of Incidence in Cross-Sectioning Application

Valery Ray^{1,2}, Oliver Zhao¹, Wen-An Chiou¹, Kirill Zaitsev³, Alexander Zaisev⁴
vray@partbeamsystech.com

⁽¹⁾ University of Maryland, College Park, MD 20742 ⁽²⁾MEO Engineering Co., Inc. Methuen MA 01844; ⁽³⁾ City College, City University of New York, NY 10031; ⁽⁴⁾ College of Staten Island, City University of New York, NY 10314

Advances in Focused Ion Beam (FIB) instrumentation and control software have transformed preparation of electron microscopy samples by Focused Ion Beam (FIB) [1] into a routine procedure, however FIB cross-sectioning remains a time-consuming process.

Dependence of etching yield on the angle of on beam incidence is known [2] and differences of etching in up-slope vs. down-slope directions were explored [3], however FIB etching at glancing-angle [4] has not yet become commonplace in practical applications.

We evaluated FIB etching with multiple-raster and single-raster scanning at glancing angle of ion beam incidence. While etching with multiple-raster scanning was not sensitive to direction of raster (*Fig. 1*), with single-raster scanning etching in down-slope direction was inhibited (*Fig. 2*), whereas up-slope direction etching by single raster scanning appeared strongly enhanced.

Cross-sectional cut (*Fig. 3*) produced by multiple etches with single-raster scanning in up-slope direction at varied angles of ion beam incidence demonstrates potential application of FIB etching at glancing angle of incidence for cross-sectioning and TEM sample preparation.

[1] R. Langford and A. Petford-Long “Preparation of transmission electron microscopy cross-section specimens using focused ion beam milling” J. Vac. Sci. Technol. A 19, 2186 (2001)

[2] Xin Xu et. al., “Focused Ion Beam induced deposition and ion milling as a function of angle of ion incidence” J. Vac. Sci. Technol. B 10, 2675 (1992)

[3] L. Gianuzzi et. al., “Optimization of Nano-Machining with Focused Ion Beams” NSTI-Nanotech Vol. 2, 2005

[4] V. Ray “FIB at Glancing Angle of Incidence” Presentation on FIB User Group at International Symposium for Testing and Failure Analysis ISTFA 2013, http://www.academia.edu/5021132/FIB_at_glancing_angle_of_incidence

FIB work performed at NISP Lab, NanoCenter was partially supported by NSF-MRSEC (DMR 05-20471, Shared Experimental Facility) and UMD

FIB work performed at College of Staten Island Nanotechnology Lab was partially supported by MEO Engineering Co., Inc. and CSI/CUNY




Figure 1: Volume of material removed by same-dose etching with multiple-raster scanning as function of direction of the raster and angle of ion beam incidence.




Figure 2: Enhancement and suppression of etching at 60 degrees angle of ion beam incidence by single-raster strategy as function of direction. All parameters are same, except of (a) 7 boxes were etched down-slope and 2 up-slope, vs. reversed direction (b) where 7 boxes were etched up-slope and 2 down-slope; (c) no tilt.




Figure 3: Enhanced cross-section by multiple single-raster up-slope cuts at varied angles of ion beam incidence (a) was superior in depth, time, and volume of removed material, as compared to single-raster cut (b) with same total ion dose at normal angle of incidence.