

Unveiling the future of precise single ion implantation

M. Hrabovsky

TESCAN GROUP, Brno, 62300, Czech Republic
milos.hrabovsky@tescan.com

J. Silvent, J. Renaud, M. Reveillard, A. Delobbe


ORSAY PHYSICS, 95 Avenue des Monts Auréliens, 13710 Fuveau, France

The Focused Ion Beam (FIB) has emerged over the past decade as an exceptionally powerful tool that has taken its place among other material science instruments. Whether configured as a single-beam or dual-beam microscope, this instrument can be used for a variety of purposes including ion microprobe, secondary ion mass spectroscopy, ion microscopy, lithography, microfabrication, ion beam etching, ion beam deposition, or ion implantation.


Combined with a Wien Filter capable of separating the different species of a Liquid Metal Alloy Ion Source (LMAIS) or plasma source, such a filtered FIB (ExB FIB) enables the use of a vast majority of ion species [1] and increases, even more, the potential applications.

ORSAY PHYSICS, part of the TESCAN GROUP, is going to introduce a new dedicated single ion implantation tool, including an *in-situ* heating stage, and that can give access to a wide range of species for the implantation application with precise control of the vertical and lateral position of the implanted ion.


[1] L. Bischoff, et al. (2016), Appl. Phys. Rev., 3 021101.



(a)



(b)



(c)

Figure 1: Photoluminescence maps acquired at a wavelength of 638 nm from two ellipses implanted at room temperature in red (a) or 800 °C in green (b) with the implantation tool. PL spectra in this spectral region are plotted in (c) to compare their intensity